

# cABR/FFR Acquisition

SPECIFICS ON RECORDING AND NOTATING


# Transducers

- Insert earphones are used to present stimuli.
- Shielded inserts are recommended.
- Personally, I have not had much luck with shielded, so I use unshielded.



# Electrodes and Placement

- Silver/Silver Chloride are recommended.
- Gold is viable alternative.
  - I used cup electrodes only at first, but now prefer clip-ons for earlobes, and cup electrodes for Cz and Ground.
- Earlobe placement is preferred for inverting
- Cz is preferred for non-inverting
- Fp1 for Ground
  - Mastoid placement is dispreferred because of PAM artifacts.



# Impedance

- Impedance should be less than 3 kOhms.
  - Best recordings when all electrodes are balanced. If necessary, accept slightly higher impedance with balance, over lower impedance on some electrodes and higher on one or two.
- I find that adherence improves with time. I may have borderline impedance while I do clicks, but it will have improved by time I am ready to measure cABR/FFR.
- For me, forehead is the hardest to keep on: tends to pull away, so I need to check it more often.



# Patient State: typically, awake



- cABR/FFR assumes an awake testee. The stimuli are typically at 80 dB SPL which means that testees usually have difficulty sleeping.
- But responses can be recorded from sleeping testee.
  - If too much artifact from muscle tension, it helps to ask testee to close eyes and try to sleep.
- Most people become “entranced” when watching a video, so movement is minimal.
  - I use an Ikea Poang chair and Poang ottoman. Most people find it very comfortable. Upside is inexpensive with low maintenance. I cover it with towels that can easily be removed and washed. Downside is that it cannot be fully reclined.
- With a very active child, it may be more time-efficient to collect fewer samples, but collect more sweeps (6 sweeps of 1024 instead of 3 sweeps of 2048).

# Settings

Montage: Cz, A1, with forehead as ground

|                      | Click                              | 40 msec DA                      | 170 msec GA and BA              |
|----------------------|------------------------------------|---------------------------------|---------------------------------|
| Gain                 | 100,000                            | 100,000                         | 100,000                         |
| Filter:<br>high pass | 100,000 Hz                         | 50 Hz                           | 50 Hz                           |
| Filter: low pass     | 1,500K                             | 3000 Hz                         | 3000 Hz                         |
| Line Filter          | on                                 | on                              | on                              |
| Rejection            | 35.0 $\mu$ V                       | 35.0 $\mu$ V                    | 35.0 $\mu$ V                    |
| Rate                 | 31.1/second                        | 10.9/sec                        | 4.35/sec                        |
| Stimulus             | 100 $\mu$ s click                  | 40 msec DA                      | 170 msec GA or BA               |
| Sweeps               | 1024 x 2                           | 2048 x 3                        | 2048 x 3                        |
| Phase                | rarefaction                        | alternating                     | alternating                     |
| Intensity            | 70 dB HL                           | 80 dB SPL                       | 80 dB SPL                       |
| Ear                  | Right and left<br>simultaneous     | Right, then Left                | Right only                      |
| Earphones            | ER3 insert earphones<br>(shielded) | ER3 insert earphones (shielded) | ER3 insert earphones (shielded) |
| Offline-filtering    | none                               | 100-1500 Hz                     | none                            |

# Worksheet: (made in WORD)

Easy to revise

Includes norms

## cABR Worksheet

|       |                |                  |                 |
|-------|----------------|------------------|-----------------|
| Name: | Birthdate/Age: | Evaluation Date: | IHS File Number |
|-------|----------------|------------------|-----------------|

### 1. Click ABR: 2 runs, 1024 each at 31.1/sec at 80dB SPL, run R & L simultaneous

|   |                           |   |
|---|---------------------------|---|
| Click Latency                           | Latency V (under 6.26 ms) | 4-7 msec: CC 1 <sup>st</sup> to 2 <sup>nd</sup> (no less than 0.89) |
| Right                                   |                           |   |
| Left                                    |                           |   |
| CC: R to L (4 to 7 msec), Cutoff = 0.89 |                           |   |

### 2. 40 msec /DA/: 3 runs of 2048 each, at 10.9/sec at 80 dB SPL, Off-lined filtered 100-1500 Hz

|          |   |   |   |   |   |   |
|----------|---|---|---|---|---|---|
| Norms:   | V | A | D | E | F | O |
| Expected |   |   |   |   |   |   |
| Right    |   |   |   |   |   |   |
| Left     |   |   |   |   |   |   |

### 3. Slope for Wave V/A of 40 msec /DA/

|                        |                    |                  |                |  |
|------------------------|--------------------|------------------|----------------|--|
| Ear                    | Wave V/A amplitude | Wave V/A latency | Wave V/A slope | (6-42 ms): CC: 1 <sup>st</sup> - 3 <sup>rd</sup> |
| Expected: greater than | At least 0.30      | Less than 1.40   | At least 0.28  | At least 0.70                                    |
| Right                  |                    |                  |                |  |
| Left                   |                    |                  |                |  |
| CC: R /DA/ to L /DA/   |                    |                  |                |  |

### 4. Pitch Tracking: 3 runs (2048) /GA/ and 3 runs (2048) each /BA/

| Syllable | Waveform # | Pitch Error    | Pitch Strength | S-R Highest Correlation | At (Latency) | CC: 1 <sup>st</sup> to 3 <sup>rd</sup> (15-60 msec) | CC: 1 <sup>st</sup> to 3 <sup>rd</sup> (60-150) |
|----------|------------|----------------|----------------|-------------------------|--------------|---|---|
| Cutoff   |            | Less than 1.73 | At least 0.43  | At least 0.43           |              | At least 0.66                                       | At least 0.55                                   |
| /GA/     |            |                |                |                         |              |   |   |
| /BA/     |            |                |                |                         |              |   |   |

### 5. Cross-phaseogram: Compare /GA/ to /BA/. Save as .jpg

10-50 msec: (bright red splotch) \_\_\_\_\_

60-170 msec: (green) \_\_\_\_\_

### 6. Save copy of waveforms as .pdf

# Demographics

Important to have the file Number because it is encrypted in the “IHS Data File”.

| cABR Worksheet |                |                  |                 |
|----------------|----------------|------------------|-----------------|
| Name:          | Birthdate/Age: | Evaluation Date: | IHS File Number |

Age is important because DA norms based on age.



# Data Collection

- Approximately 1 hour
- Collected in sound booth, but since it is a supra-threshold signal, environment can be modified: i.e. door open
- I administer behavioral tests first, then electrophysiological.
- When AEPs collected, send client home.
- Analysis done off-line, usually after client is gone.

# Click ABR

I like to include instructions on the worksheet.

Cross-correlations have to be done uniformly.


1. **Click ABR:** 2 runs, 1024 each at 31.1/sec at 80dB SPL, run R & L simultaneous

| Click Latency                           | Latency V (under 6.26 ms) | <b>4-7 msec: CC 1<sup>st</sup> to 2<sup>nd</sup> (no less than 0.89)</b> |
|---|---------------------------|--|
| Right                                   |                           |  |
| Left                                    |                           |  |
| CC: R to L (4 to 7 msec). Cutoff = 0.89 |                           |  |

# Click ABR


- Place right ear sweeps on screen
- Calculate Response Consistency (4-7 ms)
  - Correlation under 0.89 is outside of normal limits.
- Add right ear sweeps
- Mark Wave V of the summed waveform
- In over 500 clients, I have never had one without Wave V on click ABR
- Click Wave V latency and consistency is typically excellent, but there are clients with delayed responses or poor consistency.
- Repeat above steps for left ear

# Verify click ABR



# Clicks: variation

(Note: some APDs had negative correlations between first and second runs)



1. Golden Ears (blue)  
cluster tightly.

2. APDs scatter:  
Most cluster for latency, but  
scatter for correlation.

Still, abnormal latency is seen  
too.



# 40 msec DA


2. **40 msec /DA/**: 3 runs of 2048 each, at 10.9/sec at 80 dB SPL. *Off-lined filtered 100-1500 Hz.*

| Norms:   | V | A | D | E | F | O |
|----------|---|---|---|---|---|---|
| Expected |   |   |   |   |   |   |
| Right    |   |   |   |   |   |   |
| Left     |   |   |   |   |   |   |


3. Slope for Wave V/A of 40 msec /DA/

| Ear                       | Wave V/A amplitude | Wave V/A latency     | Wave V/A slope | (6-42 ms): CC:<br>1 <sup>st</sup> – 3 <sup>rd</sup> |
|---------------------------|--------------------|----------------------|----------------|---|
| Expected:<br>greater than | At least 0.30      | Less than 1.40       | At least 0.26  | At least 0.70                                       |
| Right                     |                    |                      |                |   |
| Left                      |                    |                      |                |   |
|                           |                    | CC: R /DA/ to L /DA/ |                |   |

# Measure 40 msec /DA/



# Latency of specific waveforms marked




# Peak Marking for 40 ms DA

- Activate “show cursors”
- Smart EP allows latency/amplitude markings with cursors.
  - Mainly important with Wave V/A complex.
  - Amplitude of individual waves other than V/A not measured.
- Activate all sweeps on screen.
- Under Process: Add All on Page. Move summed wave down. It is active.
- Under Process: Filter: Bandwidth. High Pass = 100; Low Pass = 1000
- Under Process: Filter Active
  - Measure latency Waves V, A, D, E, F, O. Right and left.
- Calculate slope
- Calculate Response Consistency for DA (6 to 42 ms).




# Normal 40 ms DA (2 of 3 rule)





# Response Consistency: 40 msec DA

normals: tight cluster; abnormals: wide scatter



# 170 msec GA and BA

4. Pitch Tracking: 3 runs (2048) /GA/ and 3 runs (2048) each /BA/

| Syllable | Waveform # | Pitch Error    | Pitch Strength | S-R Highest Correlation | At (Latency) | CC: 1 <sup>st</sup> to 3 <sup>rd</sup> (15-60 msec) | CC: 1 <sup>st</sup> to 3 <sup>rd</sup> (60-150) |
|----------|------------|----------------|----------------|-------------------------|--------------|---|---|
| Cutoff   |            | Less than 1.73 | At least 0.43  | At least 0.43           |              | At least 0.66                                       | At least 0.55                                   |
| /GA/     |            |                |                |                         |              |   |   |
| /BA/     |            |                |                |                         |              |   |   |

5. Cross-phaseogram: Compare /GA/ to /BA/. Save as .jpg

10-50 msec: (bright red splotch) \_\_\_\_\_

60-170 msec: (green) \_\_\_\_\_


6. Save copy of waveforms as .pdf

- Whereas the DA seems less prone to artifact, the GA and BA are more affected.

# 180 ms waveforms: GA first, then BA

- Place all sweeps on screen
- Measure Response Consistency for:
  - 1. Onset: 15-60 ms
  - 2. Vowel: 70-150 ms
- Add all sweeps, mark summed waveform (now active)
- Under “cABR” change filter settings (high pass: 75; low pass: 500). This changes Matlab module.
- Under “cABR” choose Display.
  - Find highest correlation in 6 ms to 10 ms range (absolute value)
  - Note waveform number

# Normal CC, GA 15 to 60 ms




# 180 ms GA and BA

- Under “cABR” choose Pitch-tracking
  - Find client’s file
  - Choose summed waveform for analysis
- Matlab module calculates Pitch-Tracking info and outputs a graphic
- Under “cABR” choose Crossphaseogram
  - Choose summed waveform for GA, then summed waveform for BA.
  - Matlab module calculates Cross-phaseogram for 15 to 60 ms.
  - Change Matlab module to 70-150, remeasure.



# Normal Pitch Tracking



# Pitch Tracking


- Pitch Track is plotted in the lower left of the graphic:
  - Black line represents the fundamental periodicity of the stimulus [pitch].
  - The 180 BA and 180 GA syllables have “flat” pitch that ranges from 90-110 Hz. *This is the easiest pitch to track.* There are other syllables that have distinct variation in pitch that have been used in some of the Lab’s research studies, but are not used on a regular basis.
- Yellow boxes represent the response’s periodicity. This has been time-shifted to line up with the stimulus.
- It is typical for FFR to closely follow the periodicity of the evoking stimulus, and it can be trained.

# Pitch Tracking

- The yellow boxes show how much the response's periodicity varies from the stimulus.
  - The value is given in Hz.
  - Normal Pitch Error (from my norms) is 1.72 Hz or less.
  - This value is an absolute number. It does not matter whether response was sharp or flat versus stimulus.
- Correlation Coefficient is Pearson's  $r$  between stimulus-track and response-track.
- Pitch Strength is the mean of the correlations at the thin black line. Higher is better. Normal is high.



# Cross-phaseogram




# Crossphaseogram


- Speech sound contrasts found in input stimuli are preserved in the phase of responses.
- The crossphaseogram shows phase shift as a function of time and frequency for the comparison
  - X-axis: time
  - Y-axis: frequency
  - Color-axis: Phase shift in radians



# Cross-phaseogram



IHS allows these graphics  
To be saved as .jpg.  
I add them to the report.  
Not so easy to add the  
Waveforms to report.



# Comparing GA to BA

- Expect significant phase shift in the 20-60 ms range [this is the formant transition period of the stimulus] because the syllables have different consonants and there should be a difference (color change: should be red)
- Do not expect a phase shift in the 60-170 ms range because both syllables have the same vowel (no color change: should be green).

# Completed Worksheet: helpful to highlight abnormal results

**cABR Worksheet**

|  |                                     |                                    |   |
|--|-------------------------------------|------------------------------------|---|
|  | Birthdate/Age:<br><u>12-19-2003</u> | Evaluation Date:<br><u>2-26-15</u> | IHS File Number<br><u>2015</u><br><u>2P01</u> |
|--|-------------------------------------|------------------------------------|---|

1. Click ABR: 2 runs, 1024 each at 31.1/sec at 80dB SPL, run R & L simultaneous

| Click Latency                           | Latency V (under 6.26 ms) | 4-7 msec: CC 1 <sup>st</sup> to 2 <sup>nd</sup> (no less than 0.89) |
|---|---------------------------|---|
| Right                                   | <u>5.90</u>               | <u>.822</u>   |
| Left                                    | <u>5.90</u>               | <u>.987</u>   |
| CC: R to L (4 to 7 msec). Cutoff = 0.89 |                           | <u>.974</u>   |

2. 40 msec /DA/: 3 runs of 2048 each, at 10.9/sec at 80 dB SPL, Off-lined filtered 100-1500 Hz.

| Norms:   | V           | A           | D            | E            | F            | O            |
|----------|-------------|-------------|--------------|--------------|--------------|--------------|
| Expected | <u>6.91</u> | <u>8.04</u> | <u>23.35</u> | <u>31.51</u> | <u>40.03</u> | <u>48.52</u> |
| 10 Right | <u>6.15</u> | <u>7.15</u> | <u>21.65</u> | <u>30.20</u> | <u>39.10</u> | <u>51.75</u> |
| Left     | <u>6.08</u> | <u>6.98</u> | <u>21.28</u> | <u>30.20</u> | <u>39.53</u> | <u>50.95</u> |

3. Slope for Wave V/A of 40 msec /DA/

| Ear                    | Wave V/A amplitude | Wave V/A latency | Wave V/A slope | (6-42 ms): CC: 1 <sup>st</sup> - 3 <sup>rd</sup> |
|------------------------|--------------------|------------------|----------------|--|
| Expected: greater than | At least 0.30      | Less than 1.40   | At least 0.26  | At least 0.70                                    |
| Right                  | <u>.34</u>         | <u>1.00</u>      | <u>.34</u>     | <u>.667</u>                                      |
| Left                   | <u>.32</u>         | <u>.90</u>       | <u>.35</u>     | <u>.665</u>                                      |
| CC: R /DA/ to L /DA/   |                    |                  |                | <u>.841</u>                                      |

4. Pitch Tracking: 3 runs (2048) /GA/ and 3 runs (2048) each /BA/

| Syllable # | Waveform  | Pitch Error    | Pitch Strength | S-R Highest Correlation | At (Latency) | CC: 1 <sup>st</sup> to 3 <sup>rd</sup> (15-60 msec) | CC: 1 <sup>st</sup> to 3 <sup>rd</sup> (60-150) |
|------------|-----------|----------------|----------------|-------------------------|--------------|---|---|
| Cutoff     |           | Less than 1.73 | At least 0.43  | At least 0.43           |              | At least 0.66                                       | At least 0.55                                   |
| /GA/       | <u>11</u> | <u>1.11</u>    | <u>.97</u>     | <u>.445</u>             | <u>9.00</u>  | <u>.845</u>   | <u>.904</u>                                     |
| /BA/       | <u>12</u> | <u>.72</u>     | <u>.97</u>     | <u>.489</u>             | <u>8.85</u>  | <u>.880</u>   | <u>.910</u>                                     |

5. Cross-phaseogram: Compare /GA/ to /BA/. Save as .jpg

10-50 msec: (bright red splotch) yes

60-170 msec: (green) yes

6. Save copy of waveforms as .pdf

good

1. I collected my own click norms.
2. I use the 40 msec norms from ANL (published) for latency.
3. I collected my own slope norms and cross-correlation norms.
4. Collected own norms for GA and BA.
5. I find it helpful to highlight abnormal responses on the Worksheet so I can overview easily.

# Suggested Patterns relevant to APD

- “Decoding” Deficit:
  - Poor neural encoding of transients (plosive burst) with preserved steady-state (vowel) encoding.
    - Delayed waveforms (in particular, Wave A and Wave O)
  - Poor slope for the Wave V/Wave A complex (amplitude/latency),
    - Due to either low amplitude or elongated V/A latency
- “Noise Exclusion Deficit” (Tolerance/Fading Memory):
  - Neural synchrony degraded by noise
  - Poor tracking of the fundamental frequency (Pitch Tracking)
  - Poor response consistency (ability to sustain response over time)
  - Poor stimulus-to-response correlation
- “Prosodic Deficit”:
  - Would expect FFR, Pitch-tracking, and F0 measures to be affected.